

CORDGUARD®

*The Unified System for Umbilical Cord Clamping,
Cutting & Blood Collection*

*Designed to
eliminate
spraying and
sputtering*

4-Step Actuation

*Enclosed collection reservoir
for obtaining clean,
uncontaminated blood*

*Ergonomic handle assists
in delivery of placenta*

Fully enclosed sharps

*Reliable, secure
cord clamp*

*Protective sheath keeps
tubes clean for labeling*

*Replaceable break-resistant
vacuum tubes conveniently
slide from handle when full*

A NEW OBSTETRICAL STANDARD OF CARE

Umbilical cord management that provides an excellent neonatal blood sample while minimizing exposure to blood.

OSHA
Regulations
State:

“All procedures involving blood or other potentially infectious materials shall be performed in such a manner as to minimize splashing, spraying, sputtering and generation of droplets of these substances.”

Source: OSHA 29 CFR
part 1910.1030
“Occupational
Exposure to
Bloodborne
Pathogens”

United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
Toll-free: (800) 533-4984
Tel: (801) 566-1200
Fax: (801) 566-2062

Europe

Utah Medical Products Ltd.
Athlone Business & Technology Park
Dublin Road
Athlone, County Westmeath, Ireland
Tel: (353) 90 647-3932

UTAH MEDICAL
PRODUCTS INC.



Catalog No. CRD-201
Quantity 10 per case